



macromedics

SOLUTIONS IN RADIOTHERAPY

MANUAL 126270
Instructions for use
English Instructions

MacroMedics® Respiratory Suppression Belt



MacroMedics BV

Oostbaan 670
2841ML Moordrecht
The Netherlands
+31(0)182 389 777

info@macromedics.com
www.macromedics.com



Table of Contents

1. Warnings.....	3
2. MRI Safety Information.....	3
3. Intended Use.....	3
4. Use Instructions.....	4
5. Cleaning, Storage and Disposal.....	6
6. Warranty.....	6
7. Notices.....	6
8. Glossary of Symbols.....	7

1. Warnings



- Inspect the device thoroughly before use. Do not use the device if damaged or if missing parts.
- Check individual components on MR-Safety when used together in a patient positioning system.
- The handpump and pressure gauge of the sphygmomanometer are not radiotranslucent as they contain metal parts. Therefore, keep the gauge/sphygmomanometer out of radiation pathway.
- Make sure that the device is securely fixed before positioning the patient. Place the patient carefully on the device and make sure that they lie down calmly. Instruct the patient not to reposition themselves after set-up.
- Do not overtighten the clamps, as this could result in damage to the product.
- Do not over-inflate the belt. Ensure that only the appropriate amount of pressure is applied, considering the clinical needs associated with the specific procedure being performed. Do not exceed a total pressure of 100mmHg.
- Avoid unnecessary vibrations to the gauge (e.g. falling from the table top or uncaredful storage).
- In the USA, federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.

2. MRI Safety Information



2.1 MR CONDITIONAL

The following devices are MR Conditional:

- 126270 Respiratory Suppression Belt
- 126290 Respiratory Suppression Belt HXP

A patient with this device can be safely scanned in an MR system, meeting the following conditions:

- A Static Magnetic field of 1.5 Tesla and 3 Tesla, with
- Maximum spatial field gradient of 2,500 G/cm (25T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2 W/Kg (Normal operating mode)

Due to size and location of the metallic parts, no clinically significant RF heating is expected.

3. Intended Use

Positioning of the patient during radiotherapy and radiodiagnostics, including MR where indicated.

4. Use Instructions




The **126270** Respiratory Suppression Belt can be used with all MacroMedics baseplates with an integrated rail system.




The **126290** Respiratory Suppression Belt HXP is to be used on the iBeam® couch top.

- I. Place the clamps (at either end of the belt) on the rail of the baseplate, in the desired position, and hook the clamps onto the side of the baseplate.
- II. Turn the knobs clockwise to secure the belt in place. Do not overtighten the clamps, as this could result in damage to the product.



126270 Respiratory Suppression Belt



126290 Respiratory Suppression Belt HXP


- III. Unfix the velcro connecting the top straps of the Belt and lay the top straps to the side. Position the patient in the desired position, lying on the bottom part of the belt. As part of this, instruct the patient to not lean or pull on the knobs of the device as this can damage the device.
- IV. Wrap the top straps of the belt around the patient and use the velcro to fasten the product in place. You may optionally take note of the scale (1 to 17) on the Belt to use as an indication for the (re-)placement of the strap.
- V. Ensure that the valve (Grey knob) is closed by turning it clockwise. Squeeze the handpump to inflate the belt.

Note: Do not over-inflate the belt. Ensure that only the appropriate amount of pressure is applied, considering the clinical needs associated with the specific procedure being performed. Do not exceed a total pressure of 100mmHg.

When the belt is inflated to the appropriate pressure, close the white clip with your thumb by pressing down with at least three clicks to make sure that no air can escape the tube. *



Valve



- VI. Ensure that the patient is positioned comfortably.
- VII. Place the sphygmomanometer outside of the treatment area as it contains metallic parts which are not radiotranslucent.
- VIII. To release the pressure and deflate the Respiratory Suppression Belt, first release the white clip by pushing the front end forward with your thumb. Then, turn the valve (Grey knob) counterclockwise.



- IX. To release the patient, detach the velcro straps. Instruct the patient not to pull themselves up using the knobs of the respiratory suppression belt as this can damage the device.

To release the Respiratory Suppression Belt from the baseplate, rotate the white knobs counterclockwise and carefully lift the device off the table.

***PLEASE NOTE:**

A decrease in pressure reading can be observed on the sphygmomanometer as the Belt is in use, when the white clip is closed in place. Do not worry, the pressure in the inflated cushion remains stable. What is occurring is that the part of the tube between the closed clip and the sphygmomanometer is losing pressure, but the part of the tube from the closed clip to the belt cushion remains stable at the correct pressure. This is simply because the white clip is closed in place.

After inflation and closing off of the tube with the clip, you can also detach the sphygmomanometer from the tube/belt, preventing the sphygmomanometer from damage e.g. by accidentally falling from the treatment table. To detach the sphygmomanometer from the belt, carefully pull the tube off the sphygmomanometer. During treatment it can safely be stored in its dedicated casing.



5. Cleaning, Storage and Disposal

- Clean the device before and after every use with 70% Isopropanol or soapy water, applied with a soft cloth. Rinse with water and dry properly. Follow the cleaning precautions detailed on the product.
- If unsure about the cleaning fluid, do not use it. Never use Ethanol (Ethyl alcohol), aerosol sprays, corrosive cleaning agents, solvents or abrasive detergents.
- Make sure that there is no residue of the cleaning agent left on the device prior to use.
- The devices are not suitable for sterilization or thermal disinfection.
- Having cleaned the device, perform a visual inspection of the device to ensure that it is thoroughly clean. If it is not, either repeat the cleaning process described above until the device is thoroughly clean upon visual inspection or dispose of the device to prevent the use of a soiled device.
- Store the device in a clean, ventilated environment.
- Dispose of the device in accordance with the local laws and regulations applicable to your country.

6. Warranty

The MacroMedics products are warranted to be free from defects in manufacturing, material and construction. The warranty is for 1 year starting from the date of delivery.








Product warranties remain valid provided the product was installed and/or used according to the instructions. Defects, malfunctions or failures caused by nature, misuse, abuse and unauthorised alteration or repair are not warranted.

This warranty is limited to the repair and/or replacement at the sole discretion of MacroMedics.

7. Notices

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user or patient is established.

8. Glossary of Symbols

Symbol	Symbol Name	Description of Symbol
	European Conformity	EC declaration of conformity by manufacturer.
	Manufacturer	Indicates the medical device manufacturer.
	Catalogue number	Indicates the manufacturer's catalogue number so that the medical device can be identified.
 GS1 Code: (10)	Batch Code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Consult Instruction for Use	Indicates the need for the user to consult the instructions for use.
	Caution	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.
Rx ONLY	Prescription Only	Devices for sale to physicians or trained healthcare providers only.
	MR Conditional	An item with demonstrated safety in the MR environment within defined conditions.



macromedics

SOLUTIONS IN RADIOTHERAPY

MacroMedics BV

Oostbaan 670
2841 ML Moordrecht
The Netherlands



+31 (0) 182 389 777
info@macromedics.com

www.macromedics.com

(Visit www.macromedics.com for a digital copy of this manual)



11032025V7IFU126270